

The Soft, Natural Feel Your Patients Desire,^{iv} Now with:

Xtra Projection³

More projecting than 90%
of corresponding Natrelle®
Inspira profiles by base
width.^{v,6}


Xtra Fullness³

Precision-filled with
Mentor's proprietary
cohesive MemoryGel™

Xtra Firmness³

Comparable firmness against Natrelle®
Inspira Breast Implants. †⁴

(30% higher tensile strength versus
Natrelle Inspira Breast Implant shells)^{‡5}



Women Prefer MemoryGel™ Xtra

Chosen by 9 out of 10 patients in a blinded comparison as feeling more
like natural breast than another manufacturer⁷



In a blinded comparison, surgeons chose
MemoryGel™ Xtra Breast Implants more
often than another manufacturer.⁸

v. Published projection dimensions were compared between MemoryGel™ Xtra Moderate Plus Profile vs. Inspira Moderate Profile and MemoryGel™ Xtra High Profile vs. Inspira Full Profile

iv. In-person consumer survey with 452 respondents

vi. Head-to-head blinded in-person tabletop product comparison (MemoryGel™ Xtra vs. Inspira Responsive vs. Inspira Cohesive) with 452 respondents.

head-to-head in-person tabletop product comparison (MemoryGel™ Xtra vs. Inspira Responsive vs. Inspira Cohesive) with 79 surgeons.

‡ As compared to MemoryGel® Xtra Breast Implants.

§ Head-to-head compression benchtop testing between MemoryGel® Xtra (n=4) and Natrelle Inspira (n=3). Head-to-head testing according to industry standard ASTM D412 test method for rubber properties in tension (v. 0901) between MemoryGel® (n=10) and Natrelle Inspira (n=10).

Smooth Round Breast Implants, Moderate Plus Profile Xtra



Volume	Diameter	Projection	Catalog #	Resterilizable Gel Sizer*
130 cc	8,9 cm	3,1 cm	SMPX-130	RSZ-MPX-130
160 cc	9,5 cm	3,3 cm	SMPX-160	RSZ-MPX-160
190 cc	10 cm	3,4 cm	SMPX-190	RSZ-MPX-190
215 cc	10,4 cm	3,6 cm	SMPX-215	RSZ-MPX-215
240 cc	10,8 cm	3,7 cm	SMPX-240	RSZ-MPX-240
270 cc	11,4 cm	3,8 cm	SMPX-270	RSZ-MPX-270
295 cc	11,5 cm	4 cm	SMPX-295	RSZ-MPX-295
325 cc	11,9 cm	4,1 cm	SMPX-325	RSZ-MPX-325
350 cc	12,3 cm	4,1 cm	SMPX-350	RSZ-MPX-350
370 cc	12,6 cm	4,1 cm	SMPX-370	RSZ-MPX-370
405 cc	12,7 cm	4,4 cm	SMPX-405	RSZ-MPX-405
440 cc	13,1 cm	4,5 cm	SMPX-440	RSZ-MPX-440
490 cc	13,6 cm	4,7 cm	SMPX-490	RSZ-MPX-490
545 cc	14 cm	4,9 cm	SMPX-545	RSZ-MPX-545
605 cc	14,5 cm	5,1 cm	SMPX-605	RSZ-MPX-605
645 cc	14,9 cm	5,1 cm	SMPX-645	RSZ-MPX-645
755 cc	15,7 cm	5,4 cm	SMPX-755	RSZ-MPX-755

Smooth Round Breast Implants, High Profile Xtra



Volume	Diameter	Projection	Catalog #	Resterilizable Gel Sizer*
150 cc	8,3 cm	4,1 cm	SHPX-150	RSZ-SHPX-150
180 cc	8,8 cm	4,4 cm	SHPX-180	RSZ-SHPX-180
200 cc	9,2 cm	4,5 cm	SHPX-200	RSZ-SHPX-200
235 cc	9,5 cm	4,9 cm	SHPX-235	RSZ-SHPX-235
260 cc	9,9 cm	5 cm	SHPX-260	RSZ-SHPX-260
285 cc	10,3 cm	5 cm	SHPX-285	RSZ-SHPX-285
335 cc	10,7 cm	5,4 cm	SHPX-335	RSZ-SHPX-335
355 cc	11 cm	5,5 cm	SHPX-355	RSZ-SHPX-355
380 cc	11,3 cm	5,5 cm	SHPX-380	RSZ-SHPX-380
415 cc	11,5 cm	5,8 cm	SHPX-415	RSZ-SHPX-415
450 cc	11,9 cm	5,9 cm	SHPX-450	RSZ-SHPX-450
465 cc	12,1 cm	5,8 cm	SHPX-465	RSZ-SHPX-465
490 cc	12,5 cm	5,8 cm	SHPX-490	RSZ-SHPX-490
535 cc	12,7 cm	6 cm	SHPX-535	RSZ-SHPX-535
595 cc	13,1 cm	6,3 cm	SHPX-595	RSZ-SHPX-595
650 cc	13,5 cm	6,5 cm	SHPX-650	RSZ-SHPX-650
700 cc	14,1 cm	6,5 cm	SHPX-700	RSZ-SHPX-700
755 cc	14,4 cm	6,7 cm	SHPX-755	RSZ-SHPX-755
790 cc	14,8 cm	6,7 cm	SHPX-790	RSZ-SHPX-790

SILTEX™ Round Breast Implants, Moderate Plus Profile Xtra



Volume	Diameter	Projection	Catalog #	Resterilizable Gel Sizer*
130 cc	9,1 cm	3,1 cm	TMPX-130	RSZ-MPX-130
160 cc	9,5 cm	3,3 cm	TMPX-160	RSZ-MPX-160
190 cc	10,1 cm	3,5 cm	TMPX-190	RSZ-MPX-190
215 cc	10,5 cm	3,6 cm	TMPX-215	RSZ-MPX-215
240 cc	10,9 cm	3,8 cm	TMPX-240	RSZ-MPX-240
270 cc	11,4 cm	3,9 cm	TMPX-270	RSZ-MPX-270
295 cc	11,6 cm	4,1 cm	TMPX-295	RSZ-MPX-295
325 cc	11,9 cm	4,2 cm	TMPX-325	RSZ-MPX-325
350 cc	12,2 cm	4,3 cm	TMPX-350	RSZ-MPX-350
370 cc	12,7 cm	4,3 cm	TMPX-370	RSZ-MPX-370
405 cc	12,8 cm	4,5 cm	TMPX-405	RSZ-MPX-405
440 cc	13,1 cm	4,8 cm	TMPX-440	RSZ-MPX-440
490 cc	13,6 cm	4,8 cm	TMPX-490	RSZ-MPX-490
545 cc	14,1 cm	5 cm	TMPX-545	RSZ-MPX-545
605 cc	14,4 cm	5,3 cm	TMPX-605	RSZ-MPX-605
645 cc	15 cm	5,3 cm	TMPX-645	RSZ-MPX-645
755 cc	15,7 cm	5,6 cm	TMPX-755	RSZ-MPX-755

SILTEX™ Round Breast Implants, High Profile Xtra



Volume	Diameter	Projection	Catalog #	Resterilizable Gel Sizer*
150 cc	8,4 cm	4,2 cm	THPX-150	RSZ-THPX-150
180 cc	9 cm	4,4 cm	THPX-175	RSZ-THPX-175
200 cc	9,4 cm	4,5 cm	THPX-200	RSZ-THPX-200
230 cc	9,8 cm	4,8 cm	THPX-230	RSZ-THPX-230
255 cc	10,1 cm	4,9 cm	THPX-255	RSZ-THPX-255
285 cc	10,5 cm	5,1 cm	THPX-285	RSZ-THPX-285
325 cc	10,8 cm	5,4 cm	THPX-325	RSZ-THPX-325
340 cc	11,1 cm	5,4 cm	THPX-340	RSZ-THPX-340
365 cc	11,4 cm	5,5 cm	THPX-365	RSZ-THPX-365
405 cc	11,6 cm	5,9 cm	THPX-405	RSZ-THPX-405
425 cc	11,9 cm	5,8 cm	THPX-425	RSZ-THPX-425
455 cc	12,1 cm	6 cm	THPX-455	RSZ-THPX-455
470 cc	12,3 cm	5,9 cm	THPX-470	RSZ-THPX-470
515 cc	12,6 cm	6,2 cm	THPX-515	RSZ-THPX-515
570 cc	13,2 cm	6,3 cm	THPX-570	RSZ-THPX-570
620 cc	13,4 cm	6,5 cm	THPX-620	RSZ-THPX-620
680 cc	14 cm	6,6 cm	THPX-680	RSZ-THPX-680
725 cc	14,5 cm	6,7 cm	THPX-725	RSZ-THPX-725
765 cc	14,7 cm	6,7 cm	THPX-765	RSZ-THPX-765

For further information, talk to your Mentor Sales Representative or visit mentorwllc.eu.

Information Safety Information

MENTOR® MemoryGel™ Breast Implants are indicated for breast augmentation, in women who are at least 18 years old, or for breast reconstruction. Breast implant surgery should not be performed in those women with active infection anywhere in their body, those with existing cancer or pre-cancer of their breast(s), those who have not received adequate treatment for these conditions or those who are pregnant or nursing. There are risks associated with breast implant surgery. Breast implants are not lifetime devices and breast implantation is not necessarily a one-time surgery. The most common complications with MENTOR® MemoryGel™ Breast Implants include re-operation, implant removal, capsular contracture, asymmetry, and breast pain. A lower risk of complication is implant rupture, which is most often silent. The health consequences of a ruptured silicone gel-filled breast implant have not been fully established. Screenings such as mammography, MRI, or ultrasound are recommended after initial implant surgery to assist in detecting implant rupture. Your patient needs to be informed and understand the risks and benefits of breast implants, and she should be provided with an opportunity to consult with you prior to deciding on surgery. For detailed indications, contraindications, warnings and precautions associated with the use of all MENTOR® Implantable Devices, please refer to the Product Insert Data Sheet provided with each product, or review the Important Safety Information provided at www.mentorwllc.eu. The most common complications for breast augmentation with MemoryGel™ Implants include any reoperation, capsular contracture, nipple sensation changes, and implant removal with or without replacement. The most common complications with CPG™ Implants for breast augmentation include reoperation for any reason, implant removal with or without replacement, and ptosis. A lower risk of complication is rupture. The health consequences of a ruptured silicone gel breast implant have not been fully established. MRI screenings are recommended three years after initial implant surgery and then every two years after to detect silent rupture.

1. Product Dimensions for MemoryGel and MemoryGel Xtra Breast Implants and Mentor R&D Compression Benchtop Testing - July 2017. 2. Mentor Consumer Preference Market Research Report – July 2017. 3. Product Dimensions for MemoryGel™ and MemoryGel™ Xtra Breast Implants R&D Benchtop Testing (Pich Test) - July 2017. 4. Mentor R&D Compression Benchtop Testing - July 2017. 5. Mentor Implant Shell Mechanical Properties between Allergan Natrelle and MemoryGel Xtra - August 2017. 6. Allergan Natrelle Inspira Product Catalog 2017. 7. Mentor Consumer Preference Market Research Report - July 2017. 8. Mentor Surgeon Preference Market Research Report - July 2017

All other companies' brand names mentioned herein are the trademarks of their respective owners. ©Mentor Worldwide LLC 2018 090616-180424

